

Determination of Density of Argon with a Mechanical Oscillator Densimeter

L. L. Cortés^{C, S} and E. Mendez-Lango

*División de Termometría, Laboratorio de Propiedades Termofísicas, Centro Nacional de Metrología, el
Marqués, Querétaro, México
llira@cenam.mx*

A method of calibration and characterization of an Anton Parr densimeter were investigated. Values of the density of Argon have been measured by means of a vibrating-tube densimeter at six temperatures and compared with results in the literature. The results are in good agreement. Values of the constants of calibration are obtained and densimeter is characterized.